




Medial

Innovating Healthcare, Empowering Life

Disposable Products Catalogue

www.medial.com.tr

disposet



CARDIOLOGY



We offer product range includes manifold system, manifold set, and dose control syringe, all designed to enhance precision and efficiency in cardiology. These products are essential for managing and administering medications and fluids accurately during cardiac interventions.

CARDIOLOGY

DOSE CONTROL SYRINGE

The Dose Control Syringe ensures precise delivery of the desired liquid volume, featuring finger rings for a secure grip and a luer connection for reliable attachment.

- Transparent polycarbonate body
- Accurate dosage control



MANIFOLD

Designed for angiography coronary PTCA applications, the manifold features a clear polycarbonate body capable of withstanding pressures up to 500 PSI. It includes 3-way connections and is available in both right-off and left-off variants.



MANIFOLD SET

For angiography and coronary PTCA applications, manifold set include a three-way, 500 PSI manifold with a right-angle configuration, along with extension and infusion lines, and a dose control syringe.







Innovating Healthcare Empowering Life

Address:

Bestekar Şakir Ağa Sok. Trump Villa, 6/1 Şişli - İstanbul / Turkey

E-Mail:

export@medial.com.tr

www.medial.com.tr