

I.V Cannula

- IV Cannula / IV Catheter / Branula / Intravenous Catheter
- Back-cut needle bevel designed to minimize tissue trauma
- Ported and non-porting design options
- Special rounded tip design for smooth, comfortable insertion
- High-flow performance
- Luer lock cap
- Flexible catheter resistant to kinking and bending
- Blood control feature
- International color coding for quick identification of size and specifications
- Radiopaque lines for improved visibility under imaging



Product Codes	Size	Color
14MG-IV	14G	Orange
16MG-IV	16 G	Gray
18MG-IV	18 G	Green
20MG-IV	20 G	Pink
22MG-IV	22 G	Blue
24MG-IV	24 G	Yellow

