



Suction Catheter

Standard Type

- DEHP-free and latex-free
- Smooth surface for easy insertion and enhanced glide
- Rounded, atraumatic tip
- Side eyes to help prevent tissue adhesion and blockage
- Available in different lengths and a wide range of sizes
- Color-coded funnel connector design for ease of use
- Sterile-packed
- Optimal stiffness to help prevent catheter kinking
- Length: 50 cm



Standard
Type


Product Codes	CH	Color
08MG-VC	08 Ch	Blue
10MG-VC	10 Ch	Black
12MG-VC	12 Ch	White
14MG-VC	14 Ch	Green
16MG-VC	16 Ch	Orange
18MG-VC	18 Ch	Red
20MG-SC	20Ch	Yellow



Suction Catheter

Vacuum Control Connector

- DEHP-free and latex-free
- Smooth surface for easy insertion and enhanced glide
- Rounded, atraumatic tip
- Side eyes to help prevent tissue adhesion and blockage
- Available in different lengths and a wide range of sizes
- Color-coded vacuum control connector design for ease of use and effective suction control
- Sterile-packed
- Optimal stiffness to help prevent catheter kinking
- Finger control model available
- Length: 50 cm



**Vacuum
Control
Connector**

Product Codes	CH	Color
08MG-VC	08 Ch	Blue
10MG-VC	10 Ch	Black
12MG-VC	12 Ch	White
14MG-VC	14 Ch	Green
16MG-VC	16 Ch	Orange
18MG-VC	18 Ch	Red

